





What's
new?

Europclip® PREMIUM PLUS & new AS Applicators

Ackermann®
medical

COMING
IN
2015



- ▶ No slipping or tilting of the clip in the clip applicator by incorrect leading in the magazine
Kein Verrutschen oder Kippen des Clips in der Clip Zange durch inkorrekte Führung im Magazin
- ▶ Usable with Ackermann clip applicators and appliers
Verwendbar mit Ackermann Clip Applikatoren und Zangen
- ▶ Exact guidance of the clip during clip applicator closes
Exakte Führung des Clips, während Zange schließt

AI74CS6-H

SUPER SOFT LOADING CARTRIDGE

Kleine Clips, 6 Clips/ Premium Plus Cartridge,
20 Premium Plus Cartridges/ Box, gesamt 120 Clips/ Box

*Small clips, 6 clips/ Premium Plus cartridge,
20 Premium Plus cartridges/ box, in total 120 clips/ box*

AI74CM6-H

SUPER SOFT LOADING CARTRIDGE

Mittlere Clips, 6 Clips/ Premium Plus Cartridge,
20 Premium Plus Cartridges/ Box, gesamt 120 Clips/ Box

*Medium clips, 6 clips/ Premium Plus cartridge,
20 Premium Plus cartridges/ box, in total 120 clips/ box*

AI74CML6-H

SUPER SOFT LOADING CARTRIDGE

Mittel-große Clips, 6 Clips/ Premium Plus Cartridge,
20 Premium Plus Cartridges/ Box, gesamt 120 Clips/ Box

*Medium-large clips, 6 clips/ Premium Plus cartridge,
20 Premium Plus cartridges/ box, in total 120 clips/ box*

AI74CL6-H

SUPER SOFT LOADING CARTRIDGE

Große Clips, 6 Clips/ Premium Plus Cartridge,
20 Premium Plus Cartridges/ Box, gesamt 120 Clips/ Box

*Large clips, 6 clips/ Premium Plus cartridge,
20 Premium Plus cartridges/ box, in total 120 clips/ box*

Advantages

- ▶ Resilient support also by usage with defect clip applicators
Federnder Support, auch bei Verwendung mit defekten Clip Zangen
- ▶ Loading without resistance
Widerstandlose Aufnahme
- ▶ New ergonomic design
Neues, ergonomisches Design

What's
new?

Europclip[®] PREMIUM PLUS & new AS Applicators

Ackermann[®]
medical

NEW



13-1613-105AS

Clip Applikator für kleine Europclips & Europclips Premium Plus, drehbar, doppelbeweglich, für 5,5 mm Trokare

Clip applicator for small Europclips & Europclips Premium Plus, rotatable, double action, for 5,5 mm trocars



13-1610AS

Clip Applikator für mittlere Europclips & Europclips Premium Plus, drehbar, doppelbeweglich, für 10 mm Trokare

Clip applicator for medium Europclips & Europclips Premium Plus, rotatable, double action, for 10 mm trocars



13-1615AS

Clip Applikator für mittel-große Europclips & Europclips Premium Plus, drehbar, doppelbeweglich, für 10 mm Trokare

Clip applicator for medium-large Europclips & Europclips Premium Plus, rotatable, double action, for 10 mm trocars

NEW



13-1613-106AS

Clip Applikator für große Europclips & Europclips Premium Plus, drehbar, doppelbeweglich, für 12,5 mm Trokare

Clip applicator for large Europclips & Europclips Premium Plus, rotatable, double action, for 12,5 mm trocars