

GREAT RANGE OF UNITS

CAMERAS | LIGHT SOURCES | & MORE

Cameras



FULL HD CAMERA

- ▶ real 1080p FULL HD resolution
- ▶ 16:9 widescreen
- ▶ digital zoom up to 2.5x
- ▶ automatic colour adjustment
- ▶ automatic white balance
- ▶ ...also available HD cameras, 3 chip and 1 chip cameras and more

Light Sources



LED LIGHT SOURCE

- ▶ more than 50.000 hours of life time expectancy - cost reducing
- ▶ enhanced patient safety
- ▶ brightness comparable to a 1980 W Xenon light source
- ▶ daylight colour temperature due to electronic intensity adjustment
- ▶ no dark point in the middle of the light
- ▶ ...Xenon light sources available, too

Insufflators



LAP INSUFFLATOR 45 L

- ▶ inbuilt gas pre-heating system
- ▶ external desufflation valve avoids cross-contamination
- ▶ large display shows all important measures during operation
- ▶ bottle pressure monitoring
- ▶ adjustment and pre-selection possibilities
- ▶ ...also available in 20 L and 30 L

Suction Irrigation Units



DOUBLE ROLLER PUMP

- ▶ inbuilt gas pre-heating system
- ▶ external desufflation valve avoids cross-contamination
- ▶ large display shows all important measures during operation
- ▶ bottle pressure monitoring
- ▶ adjustment and pre-selection possibilities
- ▶ ...also available in 20 L and 30 L

HF Generators



HF GENERATOR WITH HIGH-RES TOUCH SCREEN

- ▶ INCL. NEW INNOVATIVE UNIVERSAL ACKERMANN SMART DETECTION
- ▶ up to 100 memory settings
- ▶ real-time monitoring of output parameters with automatic power regulation
- ▶ voice communication possibility
- ▶ new operation modes
- ▶ ...and a lot more benefits

Shaver Systems



ARTHRO SHAVER SYSTEM

- ▶ speed adjustment from 800 - 8000 RPM
- ▶ ergonomic and light weight hand piece (autoclavable)
- ▶ silent and smooth motion
- ▶ functions commandable on control unit or footswitch
- ▶ two tools can be connected at the same time
- ▶ ...and more

... Accessories and more

OVERVIEW

- ▶ medical screens
- ▶ documentation systems
- ▶ all suitable and necessary accessories

...more laparoscopic instrumentation ins catalogue sections

16
16-2

Quality made in Germany since 1954.